


## EMG707

**Water based ferrofluid with anionic surfactant coated magnetic nano-particles**  
**MAGNETIC NANO-PARTICLES DEVELOPER KIT for Biomedical application**

Issued on Nov. 28, 2005

TN-EMG707 rev.A




### Ferrotec Corporation (JAPAN)

1-4 Midoridaira Yokaichiba  
 Chiba, 289-2131 JAPAN  
 TEL 81-479-73-6752  
 FAX 81-479-73-6602  
 ffsales@ferrotec.co.jp

### Ferrotec (USA) Corporation


40 Simon Street  
 Nashua, NH 03060 USA  
 TEL 1-603-883-9800  
 FAX 1-603-883-2308  
 sales@ferrotec.com

EMG707 is water based ferrofluid containing magnetic nano particles of iron oxide coated with anionic surfactant. The particles have a nominal diameter of about 10nm having single domain & superparamagnetic property. Therefore no hysteresis on magnetization curve can be seen as a typical data obtained by VSM (Vibrating Sampling Magnetometer). The particles also have magnetic permeability as in table and an initial susceptibility of about 0.4 typically.




VSM data for typical EMG707

The particles have about 10 nm diameter on average, however some distribution of the size can be seen as a picture of TEM (Transmittance Electron Microscope). A core of the particles are iron oxides and these are well known as compatible with living body.



XRD analysis data for typical EMG707



TEM picture for typical EMG707

### Physical properties for EMG707 (specification or typical data)

Appearance	Water based ferrofluid
Saturation magnetization of ferrofluid	9.0~11.0 mT
Average particle size (Typical data)	About 10 nm
Density of ferrofluid at 25C	1.05~1.15 g/cm <sup>3</sup>

From the X-ray analysis data by using XRD (X-Ray Diffraction spectroscopy), we can see the iron oxides are the mixture of Fe<sub>3</sub>O<sub>4</sub> and gamma-Fe<sub>2</sub>O<sub>3</sub>.

The particles can be dispersed in carrier solvent in table by just mixing. If the particles can't be dispersed completely, it can be easier by appropriate heating or ultra sound treatment.

### Typical solubility property for EMG707

water	methanol	IPA	acetone	MEK	Toluene	heptane	Xylene
OK	NG	NG	NG	NG	NG	NG	NG

Please feel free to contact Ferrotec if you need technical assistance for the particles.

### CAUTION

For Investigational Use Only. The performance characteristics of this product have not been established. This product is not sold for use in products for which prolonged contact with skin or implantation in the human body. Ferrotec Corporation and its group companies do not recommend this material as safe and effective for such uses and assumes no liability for any such use. Ferrotec Corporation and its group companies also do not warrant and specifically disclaims that the use of its ferrofluid & magnetic particles for biomedical applications covered by another party's patent do not or will not infringe the intellectual property rights of the owner of such biomedical patent.